

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

ACOT7 Antibody (Clone AT1D5)

单克隆抗体

产品货号: MA01004

第三版

**免疫原:**Recombinant human ACOT7 (1-370aa) purified from E. coli

**亚型:**IgG2a,k

**克隆:**Anti-human ACOT7 mAb, clone AT1D5, is derived from hybridization of mouse F0 myeloma cells with spleen cells from BALB/c mice immunized with a recombinant human ACOT7 protein.

**描述:**ACOT7 is a member of the acyl coenzyme family. This protein hydrolyzes the CoA thioester of palmitoyl-CoA and other long-chain fatty acids. Decreased expression of this protein may be associated with mesial temporal lobe epilepsy. Alternatively spliced transcript variants encoding distinct isoforms with different subcellular locations have been characterized.

**配方:**Liquid. In Phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol.

**浓度:**1 mg/ml

**保存:** +4 ° C 保存 (1-2 周). 长期保存在-20 ° C 或者-70 ° C. 避免反复冻融.

**应用:**ELISA, WB, ICC/IF