

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

ACTA2 Antibody (Clone AT132D3)

单克隆抗体

产品货号: MA01008

第三版

免疫原: Recombinant human ACTA2 (3-377aa) purified from E. coli

亚型: Mouse IgG1, k

克隆: Anti-human ACTA2 mAb, clone AT132D3, is derived from hybridization of mouse F0 myeloma cells with spleen cells from BALB/c mice immunized with a recombinant human ACTA2 protein.

描述: Alpha-actin-2 also known as actin, aortic smooth muscle or alpha smooth muscle actin (SMactin, alpha-SM-actin, ASMA). Actin proteins are major components of the eukaryotic cytoskeleton. At least six vertebrate actin isoforms have been identified.

配方: Liquid. In Phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol

浓度: 1 mg/ml

保存: +4 °C 保存 (1-2 周). 长期保存在 -20 °C 或者 -70 °C. 避免反复冻融.

应用: ELISA, WB, ICC/IF