

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

ADI1 Antibody (Clone AT27E8)

单克隆抗体

产品货号: MA01013

第三版

**免疫原:** Recombinant human ADI1 (1-179aa) purified from E. coli

**亚型:** IgG1,k

**克隆:** Anti-human ADI1 mAb, clone AT27E8, is derived from hybridization of mouse F0 myeloma cells with spleen cells from BALB/c mice immunized with a recombinant human ADI1 protein.

**描述:** ADI1 is an enzyme that belongs to the acireductone dioxygenase family of metal-binding enzymes, which are involved in methionine salvage. This enzyme may regulate mRNA processing in the nucleus, and may carry out different functions depending on its localization. Also, it down-regulates cell migration mediated by MMP14.

**配方:** Liquid. In Phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol.

**浓度:** 1 mg/ml

**保存:** +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.

**应用:** ELISA, WB, ICC/IF