

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Adiponectin Anti-Human Antibody (Clone 5H7)

单克隆抗体

产品货号: MA01014

第三版

免疫原: Recombinant human adiponectin (15-244aa) purified from E. coli

亚型: Mouse IgG1 heavy chain and kappa light chain

克隆: Anti-human adiponectin mAb, clone 5H7, is derived from hybridization of mouse SP2/O myeloma cells with spleen cells from BALB/c mice immunized with a recombinant human Acrp30 protein.

描述: The antibody has been tested by ELISA and Western blot analysis to assure specificity and reactivity.

配方: Liquid. In Phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol.

浓度: 1mg/ml

保存: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.

应用: ELISA, WB, ICC/IF