



本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

ADK Anti-Human Antibody (Clone AT4F8)

单克隆抗体

产品货号: MA01018

第三版

免疫原: Recombinant human ADK (22-362aa) purified from E. coli

亚型: Mouse IgG1 heavy chain and kappa light chain

克隆: AT4F8

描述: The antibody has been tested by ELISA and Western blot analysis to assure specificity and reactivity. Since application varies, however, each investigation should be titrated by the reagent to obtain optimal results. Recommended dilution range for Western blot analysis is 1:1000.

配方: Liquid. In Phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol.

浓度: 1 mg/ml

保存: +4 °C 保存 (1-2 周). 长期保存在 -20 °C 或者 -70 °C. 避免反复冻融.

应用: ELISA, WB, ICC/IF