

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

ADRM1 Antibody (Clone AT34C2)

单克隆抗体

产品货号: MA01020

第三版

**免疫原:** Recombinant human ADRM1 (1-407aa) purified from E. coli

**亚型:** IgG2b, k

**克隆:** Anti-human ADRM1 mAb, clone AT34C2, is derived from hybridization of mouse F0 myeloma cells with spleen cells from BALB/c mice immunized with a recombinant human ADRM1 protein.

**描述:** ADRM1, also known as proteasomal ubiquitin receptor ADRM1, functions as a proteasomal ubiquitin receptor. This protein is an integral plasma membrane protein which promotes cell adhesion. ADRM1 is thought to undergo O-linked glycosylation. Expression of this gene has been shown to be induced by gamma interferon in some cancer cells.

**配方:** Liquid. In Phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol.

**浓度:** 1mg/ml

**保存:** +4 ° C 保存 (1-2 周). 长期保存在-20 ° C 或者-70 ° C. 避免反复冻融.

**应用:** ELISA, WB, ICC/IF