



本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

AKR7A3 Anti-Human Antibody (Clone AT2E11)

单克隆抗体

产品货号: MA01029

第三版

免疫原: Recombinant human AKR7A3 (1-331aa) purified from E. coli

亚型: Mouse IgG1 heavy chain and kappa light chain

克隆: AT2E11

描述: Aldo-keto reductase family 7, member A 3 (AKR7A3) is a member of the aldo/keto reductase superfamily involved in the detoxification and metabolism of a variety of exogenous aldehydes and ketones.

配方: Liquid. In Phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol.

浓度: 1 mg/ml

保存: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.

应用: ELISA, WB, ICC/IF