
本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Aminoacyl tRNA synthetase complex-interacting multifunctional protein 2 (AIMP2)

单克隆抗体

产品货号: MA01075

第三版

免疫原:Denatured full length native recombinant human AIMP2.

亚型:IgG

克隆:

描述:AIMP2 is an important component for association of ARS multi-complex. AIMP2 is a potent tumor suppressor. AIMP2 inhibit cmyc activity in TGF-beta signaling, link c-IAP1 to TRAF2 by TNF-alpha stimulation and activation p53 under DNA damaged condition.

配方:Anti-AIMP2 antibody containing 0.05% sodium azide for preservative.

浓度:Please see datasheet.

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF