

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Amyloid Beta ( 11-28 ) ( 12B2 ) A Anti-Human Mouse IgG MoAb

单克隆抗体

产品货号: MA01076

第三版

免疫原:Synthetic peptide for Human Amyloid [11-28]

亚型:Mouse-Mouse hybridoma IgG1

克隆:12B2

描述:Alzheimer ' s disease (AD) is characterized by the presence of extracellular plaques and intracellular neurofibrillary tangles (NFTs) in the brain. The major protein component of these plaques is beta amyloid peptide(A<sub>β</sub>), a 40 to 43 amino acid peptide cleaved from amyloid precursor protein by beta-secretase and a putative γ-secretase.

配方:Lyophilized product from 1% BSA in PBS containing 0.05%NaN<sub>3</sub>

浓度:Please see datasheet.

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF