

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

ApoE4 (5B5) Anti-Human Mouse IgG MoAb

单克隆抗体

产品货号: MA01087

第三版

免疫原:Synthetic peptide for a part of Human ApoE4

亚型:Mouse-Mouse hybridoma

克隆:5B5

描述:Apolipoprotein-E (Apo-E) is a structural component of very low-density lipoprotein (VLDL) synthesized by the liver and intestinally synthesized chylomicrons.

配方:Lyophilized product from 1% BSA in PBS containing 0.05%NaN3

浓度:200 ug/mL

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF