

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

APP Anti-Human Antibody (Clone J4H9)

单克隆抗体

产品货号: MA01089

第三版

免疫原: Recombinant human APP (18-289aa) purified from E. coli

亚型: Mouse IgG2b heavy chain and kappa light chain

克隆: Anti-human APP mAb, clone J4H9, is derived from hybridization of mouse F0 myeloma cells with spleen cells from BALB/c mice

描述: Amyloid precursor protein (APP) is the precursor molecule whose proteolysis generates amyloid beta (A β), a 39- to 42-amino acid peptide and this amyloid fibrillar form is the primary component of amyloid plaques found in the brains of Alzheimer's disease patients.

配方: Liquid. In Phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol.

浓度: 1mg/ml

保存: +4 °C 保存 (1-2 周). 长期保存在 -20 °C 或者 -70 °C. 避免反复冻融.

应用: ELISA, WB, ICC/IF