
本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

sAPP (2B3) Anti-Human Mouse IgG MoAb

单克隆抗体

产品货号: MA01090

第三版

免疫原:Synthetic peptide in portion of C-terminus of Human sAPP (DAEFRHDSGYEVHHQK)

亚型:IgG2b

克隆:2B3

描述:Amyloid precursor protein (APP) is precursor protein of Amyloid which is major constituent of senile plaque in Alzheimer ' s disease. It is known that there are three major isoforms, APP695, APP751 and APP770, and are generated from alternative splicing of common precursor mRNA.

配方:Lyophilized product from 1% BSA in PBS containing 0.05% NaN₃

浓度:Please see datasheet.

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF