



本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Axin1 Anti-Human Antibody (At1A4)

单克隆抗体

产品货号: MA01103

第三版

免疫原: Recombinant human Axin1 (546-752aa) purified from E. coli

亚型: Mouse IgG2b heavy chain and kappa light chain

克隆: Anti-human Axin1 mAb, clone AT1A4, is derived from hybridization of mouse F0 myeloma cells with spleen cells from BALB/c mice

描述: Axin was originally identified as the product of the mouse gene called “ fused ” or fu (renamed Axin), and has since been shown to play a crucial role in controlling axis formation during embryonic development. Axin1 is characterized to be a tumor suppressor and numerous gene mutations have been identified throughout Axin1 from a number of different cancers.

配方: Liquid. In Phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol.

浓度: 1 mg/ml

保存: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.

应用: ELISA, WB, ICC/IF