

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

BAG2 Antibody (Clone AT29E9)

单克隆抗体

产品货号: MA01105

第三版

免疫原: Recombinant human BAG2 (1-211aa) purified from E. coli

亚型: IgG2a, k

克隆: Anti-human BAG2 mAb, clone AT29E9, is derived from hybridization of mouse F0 myeloma cells with spleen cells from BALB/c mice immunized with a recombinant human BAG2 protein.

描述: BCL2-associated athanogene 2, also known as BAG2, is a member of the Bag family of proteins. BAG proteins compete with Hip for binding to the Hsc70/Hsp70 ATPase domain. BAG2 is a major component of the HSC 70/CHIP chaperone-dependent ubiquitin ligase complex and acts to disrupt CHIP-mediated ubiquitylation.

配方: Liquid. In Phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol.

浓度: 1 mg/ml

保存: +4 °C 保存 (1-2 周). 长期保存在 -20 °C 或者 -70 °C. 避免反复冻融.

应用: ELISA, WB, ICC/IF