

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Bid Anti-Human Antibody (Clone 4D3)

单克隆抗体

产品货号: MA01115

第三版

**免疫原:**Recombinant human BID (1-195aa) purified from E. coli

**亚型:**Mouse IgG1 heavy chain and kappa light chain

**克隆:**Anti-human BID mAb, clone 4D3, is derived from hybridization of mouse F0 myeloma cells with spleen cells from BALB/c mice

**描述:**BID is a pro-apoptotic Bcl-2 protein containing only the BH3 domain. In response to apoptotic signaling, BID interacts with another Bcl-2 family protein, Bax, leading to the insertion of Bax into the outer mitochondrial membrane.

**配方:**Liquid. In Phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol.

**浓度:**1 mg/ml

**保存:** +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

**应用:**ELISA, WB, ICC/IF