

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

BRCC36 Anti-Human Antibody (Clone AT3B1)

单克隆抗体

产品货号: MA01122

第三版

免疫原:Recombinant BRCC36 (1-316aa) purified from E. coli

亚型:Mouse IgG1 heavy chain and kappa light chain

克隆:Anti-human BRCC36 mAb, clone AT3B1, is derived from hybridization of mouse F0 myeloma cells with spleen cells from BALB/c mice

描述:Description:..BRCC36 (BRCA1/BRCA2-containing complex, subunit 36) is a holoenzyme complex containing BRCA1, BRCA2. BRCC36 is an ubiquitin E3 ligase complex that has a role in DNA repair and enhances cellular survival following DNA damage. This protein is also thought to be involved in the cellular response to ionizing radiation and progression through the G2/M checkpoint.

配方:Liquid. In Phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol.

浓度:1 mg/ml

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF