



本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

c-kit (12A8) Anti-Human Mouse IgG MoAb

单克隆抗体

产品货号: MA01126

第三版

**免疫原:** Recombinant protein of extracellular domain of human c-Kit

**亚型:** IgG1

**克隆:** 12A8

**描述:** This monoclonal antibody recognizes c-Kit by W.B. and I.P. applications and it enables FACS analysis of c-Kit-expressing cells. It has been confirmed that the addition of this antibody does not induce ligand-binding-dependent phosphorylation of c-Kit and it may be used as a neutralizing antibody.

**配方:** Lyophilized product from 1% BSA in PBS containing 0.05%NaN<sub>3</sub>

**浓度:** 100 ug/mL

**保存:** +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

**应用:** ELISA, WB, ICC/IF