

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

CAPG Antibody (Clone AT1D10)

单克隆抗体

产品货号: MA01137

第三版

免疫原: Recombinant human CAPG (1-348aa) purified from E. coli

亚型: Mouse IgG2b, k

克隆: Anti-human CAPG mAb, clone AT1D10, is derived from hybridization of mouse F0 myeloma cells with spleen cells from BALB/c mice immunized with a recombinant human CAPG protein.

描述: Macrophage-capping protein, also known as CAPG, belongs to the gelsolin/villin family of actin-regulatory proteins. CAPG is a calcium-sensitive DNA-binding protein that plays a role in macrophage function. It is expressed in macrophages and macrophage-like cells and can localize both to the nucleus and the cytoplasm. This protein reversibly blocks the barbed ends of F-actin filaments in a Ca²⁺ and phosphoinositide-regulated manner, but does not sever preformed actin filaments.

配方: Liquid. In Phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol.

浓度: 1mg/ml

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C 或者-70 ° C. 避免反复冻融.

应用: ELISA, WB, ICC/IF