

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Cathepsin L Anti-Human Antibody (Clone 2H7)

单克隆抗体

产品货号: MA01142

第三版

免疫原:Recombinant human Cathepsin L (114-333aa) purified from E. coli

亚型:Mouse IgG2b heavy chain and kappa light chain

克隆:Anti-human Cathepsin L mAb, clone 2H7, is derived from hybridization of mouse F0 myeloma cells with spleen cells from BALB/c mice

描述:Cathepsin L, a lysosomal endopeptidase expressed in most eukaryotic cells, is a member of the papain-like family of cysteine proteinases. Cathepsin L can be induced by many different signaling events, including growth factor, secondary messenger and tumor promoters.

配方:Liquid. In Phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol.

浓度:1mg/ml

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF