



本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Ccr6 Anti-Human Antibody (Clone 4C6)

单克隆抗体

产品货号: MA01147

第三版

**免疫原:**Recombinant human CCR6 (1-46aa) purified from E. coli

**亚型:**Mouse IgG1 heavy chain and kappa light chain

**克隆:**Anti-human CCR6 mAb, clone 4C6, is derived from hybridization of mouse F0 myeloma cells with spleen cells from BALB/c mice

**描述:**CCR6 is a member of the beta chemokine receptor family, which is predicted to be a seven transmembrane protein similar to G protein-coupled receptors. The gene is preferentially expressed by immature dendritic cells and memory T cells.

**配方:**Liquid. In Phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol.

**浓度:**1mg/ml

**保存:** +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.

**应用:**ELISA, WB, ICC/IF