

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Anti-CD252/OX40L (TAG34) Mouse IgG MoAb

单克隆抗体

产品货号: MA01156

第三版

免疫原:Lysate of HTLV-1 persistently-infected cell line (MT-2)

亚型:IgG1

克隆:TAG34

描述:CD252/OX40L is a molecule identical to gp34 expressed in T cells infected with HTLV-1 and it is expressed on activated B cell, dendritic cells and vascular endothelial cells as a trimer belonging to TNF superfamily. CD252/OX40L have a role in proliferation of T cells and production of cytokines through OX40 expressed on activated T cells. In addition, OX40L expressed on the activated vascular endothelial cells are involved in activated T cell invasion to inflamed area.

配方:Lyophilized product from PBS containing 1 % BSA and 0.05 % NaN₃

浓度:Please see datasheet.

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF