

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Cellubrevin Anti-Human Antibody (clone AT4G9)

单克隆抗体

产品货号: MA01172

第三版

免疫原: Recombinant human cellubrevin (1-77aa) purified from E. coli

亚型: Mouse IgG2b heavy chain and kappa light chain

克隆: Anti-human VAMP3 mAb, clone AT4G9, is derived from hybridization of mouse F0 myeloma cells with spleen cells from BALB/c mice

描述: Cellubrevin, also known as VAMP3, is a member of the Vamp/synaptobrevin family. Vamps, syntaxins, and the 25-kD synaptosomal-associated protein are the main components of a protein complex involved in the docking and/or fusion of synaptic vesicles with the presynaptic membrane.

配方: Liquid. In Phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol.

浓度: 1mg/ml

保存: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.

应用: ELISA, WB, ICC/IF