

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

CIAO1 Antibody (Clone AT6C9)

单克隆抗体

产品货号: MA01176

第三版

免疫原: Recombinant human CIAO1 (1-339aa) purified from E. coli

亚型: IgG1,k

克隆: Anti-human CIAO1 mAb, clone AT6C9, is derived from hybridization of mouse F0 myeloma cells with spleen cells from BALB/c mice immunized with a recombinant human CIAO1 protein.

描述: CIAO1 is a essential component of the cytosolic iron-sulfur (Fe/S) protein assembly machinery. It is a multiprotein complex that mediates the incorporation of iron-sulfur cluster into extramitochondrial Fe/S proteins. CIAO1 seems to specifically modulate the transactivation activity of WT1 and may function to regulate the physiological functions of WT1 in cell growth and differentiation. As part of the mitotic spindle-associated MMXD complex, it may play a role in chromosome segregation.

配方: Liquid. In Phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol.

浓度: 1 mg/ml

保存: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.

应用: ELISA, WB, ICC/IF