



本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Clusterin Anti-Human Antibody (Clone 1A11)

单克隆抗体

产品货号: MA01184

第三版

免疫原:Recombinant human Clusterin (1-333aa) purified from E. coli

亚型:Mouse IgG1 heavy chain and kappa light chain

克隆:Anti-human Clusterin mAb, clone 1A11, is derived from hybridization of mouse SP2/O myeloma cells with spleen cells from BALB/c mice

描述:Clusterin has been implicated in numerous processes including active cell death, immune regulation, cell adhesion and morphological transformation. Clusterin (previously named SP-40) is a glycoprotein, 80kDa comprising two distinct 40kDa polypeptide chains.

配方:Liquid. In Phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol.

浓度:1mg/ml

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF