

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Cofilin 1 Anti-Human Antibody (Clone At1C1)

单克隆抗体

产品货号: MA01188

第三版

免疫原: Recombinant human Cofilin 1 (1-166aa) purified from E. coli

亚型: Mouse IgG1 heavy chain and kappa light chain

克隆: Anti-human Cofilin 1 mAb, clone AT1C1, is derived from hybridization of mouse F0 myeloma cells with spleen cells from BALB/c mice

描述: Cofilin is a family of actin-binding proteins which disassembles actin filaments. It is a widely distributed intracellular actin-modulating protein that binds and depolymerizes filamentous F-actin and inhibits the polymerization of monomeric G-actin in a pH-dependent manner. It is involved in the translocation of actin-cofilin complex from cytoplasm to nucleus.

配方: Liquid. In Phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol.

浓度: 1mg/ml

保存: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.

应用: ELISA, WB, ICC/IF