

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

COTL1 Antibody (Clone AT1D6)

单克隆抗体

产品货号: MA01190

第三版

免疫原: Recombinant human FCGR3A (18-208aa) purified from E. coli

亚型: IgG2a ,k

克隆: AT1D6

描述: FCGR3A (Low affinity immunoglobulin gamma Fc region receptor III-A), also known as CD16A, is 254 amino acid protein. This protein is involved in the removal of antigen-antibody complexes from the circulation, as well as other antibody-dependent responses. It mediates antibody-dependent cellular cytotoxicity (ADCC) and other antibody-dependent responses, such as phagocytosis. Mutations in this gene have been linked to susceptibility to recurrent viral infections, susceptibility to systemic lupus erythematosus, and alloimmune neonatal neutropenia.

配方: Liquid. In Phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol.

浓度: 1mg/ml

保存: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.

应用: ELISA, WB, ICC/IF