

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

CRABP1 Anti-Human Antibody (Clone AT1A1)

单克隆抗体

产品货号: MA01194

第三版

免疫原:Recombinant human CRABP1 (1-137aa) purified from E. coli

亚型:Mouse IgG2b heavy chain and kappa light chain

克隆:Anti-human CRABP1 mAb, clone AT1A1, is derived from hybridization of mouse F0 myeloma cells with spleen cells from BALB/c mice

描述:CRABP1 (cellular retinoic acid-binding protein 1) is a member of specific carrier proteins for members of the vitamin A family. CRABP1 is assumed to play an important role in retinoic acid-mediated differentiation and proliferation processes.

配方:Liquid. In Phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol.

浓度:1 mg/mL

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF