

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Crabp2 Anti-Human Antibody (Clone At2E11)

单克隆抗体

产品货号: MA01195

第三版

免疫原:Recombinant human CRABP2 (1-138aa) purified from E. coli

亚型:Mouse IgG2a heavy chain and kappa light chain

克隆:Anti-human CRABP2 mAb, clone AT2E11, is derived from hybridization of mouse F0 myeloma cells with spleen cells from BALB/c mice

描述:The cellular retinoic acid-binding protein II (CRABP-II) is involved in the conversion of vitamin A into its intracellular active form retinoic acid, which regulate the genes responsible for lipid metabolism and adipocyte differentiation.

配方:Liquid. In Phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol.

浓度:1 mg/mL

保存: +4 °C 保存 (1-2 周). 长期保存在-20 °C 或者-70 °C. 避免反复冻融.

应用:ELISA, WB, ICC/IF