



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

CREG1 Monoclonal Antibody (Clone AT1C6)

单克隆抗体

产品货号: MA01197

第三版

**免疫原:** Recombinant human CREG1 (32-220aa) purified from E. coli

**亚型:** IgG2b,k

**克隆:** Anti-human CREG1 mAb, clone AT1C6, is derived from hybridization of mouse F0 myeloma cells with spleen cells from BALB/c mice immunized with a recombinant human CREG1 protein.

**描述:** CREG1 activity may contribute to the transcriptional control of cell growth and differentiation. Antagonizes transcriptional activation and cellular transformation by the adenovirus E1A protein. CREG1 is broadly expressed in adult tissues and is regulated during embryonic development. CREG1 is a secreted glycoprotein which enhances differentiation of mouse embryonic stem cells and human NTERA-2 cells.

**配方:** Liquid. In Phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol.

**浓度:** 1 mg/mL

**保存:** +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.

**应用:** ELISA, WB, ICC/IF