



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

-Crystallin A Anti-Human Antibody (Clone c9F2)

单克隆抗体

产品货号: MA01205

第三版

免疫原:Recombinant human Crystallin alpha A (1-173aa) purified from E. coli

亚型:Mouse IgG1 heavy chains and light chain

克隆:Anti-human Crystallin alpha A mAb, clone c9F2, is derived from hybridization of mouse SP2/O myeloma cells with spleen cells from BALB/c mice

描述:crystallins alpha are composed of two gene products; alpha A and alpha B, for acidic and basic, respectively. crystallins alpha can be induced by heat shock and are members of the small heat shock protein (sHSP also known as the HSP20).

配方:Liquid. In Phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol.

浓度:1 mg/mL

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF