
本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

-Crystallin B Anti-Human Antibody (Clone 2E8)

单克隆抗体

产品货号: MA01206

第三版

免疫原:Recombinant human Crystallin alpha B (1-175 aa) purified from E. coli

亚型:Mouse IgG2b heavy chains and light chain

克隆:Anti-human Crystallin alpha B mAb, clone 2E8, is derived from hybridization of mouse SP2/O myeloma cells with spleen cells from BALB/c mice

描述:Crystallins alpha are composed of two gene products; Crystallin alpha A and B, for acidic and basic, respectively. Crystallins alpha can be induced by heat shock and are members of the small heat shock protein (sHSP also known as the HSP20).

配方:Liquid. In Phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol.

浓度:1 mg/mL

保存: +4 °C 保存 (1-2 周). 长期保存在-20 °C或者-70 °C. 避免反复冻融.

应用:ELISA, WB, ICC/IF