

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

CSK antibody (clone AT7C6)

单克隆抗体

产品货号: MA01207

第三版

免疫原: Recombinant human CSK (1-450aa) purified from E. coli

亚型: IgG2b, k

克隆: Anti-human CSK mAb, clone AT7C6, is derived from hybridization of mouse F0 myeloma cells with spleen cells from BALB/c mice immunized with a recombinant human CSK protein.

描述: CSK is a tyrosine-protein kinase CSK that plays an important role in the regulation of cell growth, differentiation, migration and immune response. Suppresses signaling by various surface receptors, including T-cell receptor (TCR) and B-cell receptor (BCR) by phosphorylating and maintaining inactive several positive effectors such as FYN or LCK.

配方: Liquid. In Phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol.

浓度: 1 mg/ml

保存: +4 °C 保存 (1-2 周). 长期保存在 -20 °C 或者 -70 °C. 避免反复冻融.

应用: ELISA, WB, ICC/IF