



本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

CSTA Antibody (Clone AT1B12)

单克隆抗体

产品货号: MA01208

第三版

免疫原: Recombinant human CSTA (1-98aa) purified from E. coli

亚型: IgG1, k

克隆: Anti-human CSTA mAb, clone AT1B12, is derived from hybridization of mouse F0 myeloma cells with spleen cells from BALB/c mice immunized with a recombinant human CSTA protein.

描述: CSTA, also known as Stefin A, is thiol protease inhibitors that form complexes with papain and the cathepsins B, H and L. The protein is one of the precursor proteins of the cornified cell envelope in keratinocytes and plays a role in epidermal development and maintenance. CSTA is thought to be a tumor suppressor, since an inverse correlation seems to exist between the level of CSTA and tumor progression, in several malignancies.

配方: Liquid. In Phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol.

浓度: 1mg/ml

保存: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.

应用: ELISA, WB, ICC/IF