
本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

CTBP1 Anti-Human Antibody (Clone AT4D6)

单克隆抗体

产品货号: MA01210

第三版

免疫原: Recombinant human CTBP1 (1-440aa) purified from E. coli

亚型: Mouse IgG2b heavy chain and light chain

克隆: AT4D6

描述: The C-terminal-binding protein (CtBP) was originally identified as a cellular protein that interacts with a C-terminal motif of adenovirus E1A oncoproteins.

配方: Liquid. In Phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol.

浓度: 1 mg/ml

保存: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.

应用: ELISA, WB, ICC/IF