

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

CYR61 Anti-Human Antibody (Clone3H3)

单克隆抗体

产品货号: MA01216

第三版

免疫原: Recombinant human CYR61 (168-242aa) purified from E. coli

亚型: Mouse IgG1 heavy chain and kappa light chain

克隆: Anti-human CYR61 mAb, clone 3H3, is derived from hybridization of mouse FO myeloma cells with spleen cells from BALB/c mice

描述: Cysteine-rich angiogenic inducer 61, also known as CYR61, is a growth factor-inducible, immediate-early gene that has multifaceted activities in various cancers. CYR61 is a secreted heparin-binding extracellular matrix-associated protein whose activities include the promotion of adhesion and chemotaxis, and the stimulation of fibroblast and endothelial cell growth.

配方: Liquid. In Phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol.

浓度: 1 mg/mL

保存: +4 °C 保存 (1-2 周). 长期保存在 -20 °C 或者 -70 °C. 避免反复冻融.

应用: ELISA, WB, ICC/IF