

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

FABP1 Anti-Human Antibody (Clone 2G4)

单克隆抗体

产品货号: MA01244

第三版

免疫原: Recombinant human FABP1 (1-127aa) purified from E. coli

亚型: Mouse IgG1 heavy chain and kappa light chain

克隆: Anti-human FABP1 mAb, clone 2G4, is derived from hybridization of mouse F0 myeloma cells with spleen cells from BALB/c mice

描述: FABP1 (Fatty acid binding protein1) encodes the fatty acid binding protein found in liver. FABP1 is composed of ten antiparallel strands that form a barrel with a bigger binding pocket than the other FABPs allowing it to accommodate two fatty acid.

配方: Liquid. In Phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol.

浓度: 1 mg/mL

保存: +4 °C 保存 (1-2 周). 长期保存在-20 °C 或者-70 °C. 避免反复冻融.

应用: ELISA, WB, ICC/IF