



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

FAS (14G5) Anti-Human Mouse IgG MoAb

单克隆抗体

产品货号: MA01249

第三版

**免疫原:** Synthetic peptide of human FAS

**亚型:** IgG2b

**克隆:** 14G5

**描述:** Animal fatty acid synthase (FAS) is a homodimer protein, which synthesizes long-chain fatty acids and is rich in liver, brain, breast and lung. However, the precise cellular localization of FAS in human tissue has not been elucidated

**配方:** Lyophilized product from 1% BSA in PBS containing 0.05%NaN<sub>3</sub>

**浓度:** 100 ug/ml

**保存:** +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.

**应用:** ELISA, WB, ICC/IF