



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Galectin-3 (38B2) Anti-Human Mouse IgG MoAb

单克隆抗体

产品货号: MA01266

第三版

免疫原:Recombinant Galectin-3 protein

亚型:IgG1

克隆:38B2

描述:Galectin is widely distributed in nematodes, insects, and porifers, as well as vertebrates, and it has also been found to be present in true fungi. Galectin does not just occur in the cytoplasm, it is also present in the nucleus, on the cell surface, in the extracellular matrix, etc., and it is thought to be involved in many biological phenomena, including development, differentiation, morphogenesis, tumor metastasis, cell death, and RNA splicing.

配方:Lyophilized product from 1% BSA in PBS containing 0.05% NaN<sub>3</sub>

浓度:100 ug/mL

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF