



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

GAS7 Anti-Human Antibody (Clone AT4H8)

单克隆抗体

产品货号: MA01273

第三版

免疫原:Recombinant human GAS7 (1-416aa) purified from E. coli

亚型:Mouse IgG1 heavy chain and kappa light chain

克隆:Anti-human GAS7 mAb, clone AT4H8, is derived from hybridization of mouse F0 myeloma cells with spleen cells from BALB/c mice

描述:Growth arrest-specific 7, also known as GAS7, is expressed primarily in terminally differentiated brain cells and predominantly in mature cerebellar Purkinje neurons. GAS7 plays a putative role in neuronal development by promoting maturation and morphological differentiation of cerebellar neurons.

配方:Liquid. In Phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol.

浓度:1 mg/ml

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF