

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

Geminin Anti-Human Antibody (Clone 2H7)

单克隆抗体

产品货号: MA01277

第三版

免疫原: Recombinant human Geminin (1-209 aa) purified from E. coli

亚型: Mouse IgG1 heavy chain and kappa light chain

克隆: Anti-human Geminin mAb, clone 2H7, is derived from hybridization of mouse SP2/0 myeloma cells with spleen cells from BALB/c mice

描述: Geminin is a 25 kDa protein, which inhibits DNA replication and is degraded during the mitotic phase of the cell cycle. Geminin controls replication by binding to the licensing factor Cdt1, and is involved in neural differentiation.

配方: Liquid. In Phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol.

浓度: 1 mg/ml

保存: +4 °C 保存 (1-2 周). 长期保存在 -20 °C 或者 -70 °C. 避免反复冻融.

应用: ELISA, WB, ICC/IF