



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

GFAP (G-25-8-3) Anti-Bovine Mouse IgG MoAb

单克隆抗体

产品货号: MA01278

第三版

免疫原:Bovine GFAP

亚型:IgG2b

克隆:G-25-8-3

描述:GFAP (Glial fibrillary acidic protein) is acidic protein which has molecular weight of 50kDa, and the protein is the principal constituent of glial fibrillary acidic protein. As GFAP exists in astroglia specifically, this antibody can be used as marker which recognize astroglia specifically.

配方:Lyophilized product containing 0.05%NaN<sub>3</sub>

浓度:Please see datasheet.

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF