



本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

GPC4 Antibody (Clone AT51E3)

单克隆抗体

产品货号: MA01287

第三版

免疫原: Recombinant human GPC4 (401-529aa) purified from E. coli

亚型: Mouse IgG2a, k

克隆: AT51E3

描述: The Glypican (GPC) family of cell surface heparan sulfate proteoglycans is expressed in a tissue specific and developmentally regulated fashion. Expression of Glypican 4 has been reported in the mouse brain at embryonic day 10 and later stages and low during early tubule formation but is up-regulated in mature tubules. Accordingly it is not detected in the midline of the embryo at the initial stage of neural tube closure, suggesting that GPC4 is unlikely to play an essential role in convergent extension in the mouse.

配方: Liquid. In Phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol.

浓度: 1 mg/ml

保存: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.

应用: ELISA, WB, ICC/IF