

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

GPS2 Anti-Human Antibody (Clone 1E2)

单克隆抗体

产品货号: MA01290

第三版

免疫原: Recombinant human GPS2 (161-311aa) purified from E. coli

亚型: Mouse IgG3 heavy chain and kappa light chain

克隆: Anti-human GPS2 mAb, clone 1E2, is derived from hybridization of mouse SP2/0 myeloma cells with spleen cells from BALB/c mice

描述: GPS2(G-protein pathway suppressor 2), also called AMF1, is a human nuclear protein of 327 amino acids involved in G protein-mitogen-activated protein kinase(MAPK) signaling cascades.

配方: Liquid. In Phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol.

浓度: 1 mg/ml

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用: ELISA, WB, ICC/IF