

---

本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

---

Growth Hormone Anti-Human Antibody (Clone g3H5)

单克隆抗体

产品货号: MA01294

---

第三版

**免疫原:** Recombinant human growth hormone (27-217aa) purified from E. coli

**亚型:** Mouse IgG2b heavy chain and kappa light chain

**克隆:** Anti-human growth hormone mAb, clone g3H5, is derived from hybridization of mouse SP2/O myeloma cells with spleen cells from BALB/c mice

**描述:** Human growth hormone (hGH, somatotropin), a protein hormone from the anterior lobe of the pituitary gland, comprises a single polypeptide chain of 191 amino acid residues.

**配方:** Liquid. In Phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol.

**浓度:** 1 mg/ml

**保存:** +4 °C 保存 (1-2 周). 长期保存在 -20 °C 或者 -70 °C. 避免反复冻融.

**应用:** ELISA, WB, ICC/IF