



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

GSK3 Anti-Human Antibody (Clone 1F7)

单克隆抗体

产品货号: MA01296

第三版

免疫原: Recombinant human GSK3 (341-420aa) purified from E. coli

亚型: Mouse IgG1 heavy chain and kappa light chain

克隆: Anti-human GSK3 mAb, clone 1F7, is derived from hybridization of mouse SP2/O myeloma cells with spleen cells from BALB/c mice immunized

描述: Glycogen synthase kinase 3 (GSK3) is a serine/threonine kinase involved in insulin, growth factor and Wnt signaling. Glycogen synthase kinase 3 (GSK3) is a molecular regulator of cell fate in the central nervous system and a target of the phosphatidylinositol 3-kinase (PI3K) pathway.

配方: Liquid. In Phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol.

浓度: 1 mg/ml

保存: +4 °C 保存 (1-2 周). 长期保存在 -20 °C 或者 -70 °C. 避免反复冻融.

应用: ELISA, WB, ICC/IF