

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

Anti-GST Antibody (Clone 1E5)

单克隆抗体

产品货号: MA01299

第三版

**免疫原:** Recombinant Glutathione S transferase (GST) purified from E. coli

**亚型:** Mouse IgG2b heavy chain and kappa light chain

**克隆:** Anti-GST mAb, clone 1E5, is derived from hybridization of mouse SP2/0 myeloma cells with spleen cells from BALB/c mice

**描述:** Glutathione S transferases (GSTs) are a family of enzymes that play an important role in detoxification by catalyzing the conjugation of many hydrophobic and electrophilic compounds with reduced glutathione. GST proteins can be purified by immobilized glutathione affinity chromatography.

**配方:** Liquid. In Phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol.

**浓度:** 1 mg/ml

**保存:** +4 °C 保存 (1-2 周). 长期保存在 -20 °C 或者 -70 °C. 避免反复冻融.

**应用:** ELISA, WB, ICC/IF