
本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

H1N1/HA1 Anti-Human Antibody (Clone AT1G7)

单克隆抗体

产品货号: MA01302

第三版

免疫原: Recombinant human H1N1/HA1 (18-344aa) purified from Baculovirus

亚型: Mouse IgG1 heavy chain and light chain

克隆: 1G7

描述: Influenza virus, an enveloped virus of the Orthomyxoviridae family, has a unique capacity for genetic variation that is based in two molecular features of the virus family.

配方: Liquid. In Phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol.

浓度: 1mg/ml

保存: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.

应用: ELISA, WB, ICC/IF