

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

H5N1/HA1 Anti-Human Antibody (Clone AT2B7)

单克隆抗体

产品货号: MA01304

第三版

免疫原:Recombinant human H5N1/HA1 (17-338aa) purified from Baculovirus

亚型:Mouse IgG1heavy chain and kappa light chain

克隆:2B7

描述:H5N1 is a subtype of the species Influenza A virus of the Influenzavirus A genus of the Orthomyxoviridae family. consists of single-stranded eight-segment negative-sense genomic RNAs, helical viral ribonucleoprotein (RNP) complexes (RNA segments NP, PB2, PB1 and PA), three viral envelope proteins (hemagglutinin [HA], neuraminidase [NA], and M2 ion channel), and a matrix (M1) protein.

配方:Liquid. In Phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol.

浓度:1mg/ml

保存: +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

应用:ELISA, WB, ICC/IF