

本公司提供的电子版本说明书仅供参考，实验请以收到的纸质手册为准。

STAT1 Antibody (Clone AT1F4)

单克隆抗体

产品货号: MA01312

第三版

**免疫原:** Recombinant human STAT1 (1-712aa) purified from E. coli

**亚型:** IgG1,k

**克隆:** Anti-human STAT1 mAb, clone AT1F4, is derived from hybridization of mouse F0 myeloma cells with spleen cells from BALB/c mice immunized with a recombinant human STAT1 protein.

**描述:** STAT1, also known as STAT91, is a member of the STAT protein family. In response to cytokines and growth factors, STAT family members are phosphorylated by the receptor associated kinases, and then form homo- or heterodimers that translocate to the cell nucleus where they act as transcription activators. STAT1 is activated by peptide hormones, growth factors and cytokines, particularly cytokines of the interferon family.

**配方:** Liquid. In Phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol.

**浓度:** 1 mg/ml

**保存:** +4 ° C 保存 (1-2 周). 长期保存在-20 ° C或者-70 ° C. 避免反复冻融.

**应用:** ELISA, WB, ICC/IF