



本公司提供的电子版说明书仅供参考，实验请以收到的纸质手册为准。

VHL Anti-Human Antibody (Clone 52A11)

单克隆抗体

产品货号: MA01371

第三版

免疫原: Recombinant human VHL (1-154aa) purified from E. coli

亚型: Mouse IgG2b heavy chain and kappa light chain

克隆: 52A11

描述: Von Hippel-Lindau disease (VHL) is a dominant inherited syndrome characterized by the predisposition to develop various kinds of benign and malignant tumors, including clear cell renal carcinomas, pheochromocytomas and hemangioblastomas of the central nervous system and retina.

配方: Liquid. In Phosphate-Buffered Saline (pH 7.4) with 0.02% Sodium Azide, 10% Glycerol.

浓度: 1mg/ml

保存: +4 ° C 保存 (1-2 周). 长期保存在 -20 ° C 或者 -70 ° C. 避免反复冻融.

应用: ELISA, WB, ICC/IF